

Scan Bodies, Digital Libraries and Titanium Bases



Guidelines to Download and Install the PREAT Library Files for 3Shape Titanium Bases

Step 1. Download files for 3Shape at
www.preat.com/portfolio/digital-libraries/



3Shape


Digital Library Download for 3Shape



3Shape Guideline to Download and Install



Step 2. Save file to a location on the computer that is easily accessible.



Step 3. Open 3Shape Dental System Control Panel and select Import/Export in the Tools menu.



Dental System Control Panel

Home Back Save Help

Home Page

- System Settings**
 - System Settings
 - Design options
 - Services
- DentalManager**
 - General
 - Manufacturing Inbox
 - ERP and processing time
 - Order import
 - Working days
 - 3rd party applications
 - Milling machines
- Site Settings**
 - This site
 - Dentists
 - Labs
 - Manufacturers
 - Operators
 - Countries
 - Manufacturing processes
- Tools**
 - Subscription Management (Dongle)
 - Import/Export**
 - Download Center
 - Training Center
 - 3Shape Communicate

- Basic elements**
 - Materials
 - Colors
 - Margin line
 - Attachments
 - 2D Design overlays
 - CAD blocks
- Anatomy elements**
 - ScanIt library
 - Anatomy and Pontic libraries
 - Smile libraries
 - Artificial teeth
 - Crowns
 - Crown Pontics
 - Inlays
 - Onlays/Veneers
 - Temporary Crowns
 - Temporary Pontics
 - Temporaries on prepared model
 - Tabletops
- Frame elements**
 - Copings
 - Frame Pontics
 - Waxups
 - Primary telescopes
- Abutments**
 - Top cap libraries
 - Implant systems

- Post and Core**
 - Post and Core systems
 - Post and Cores
- Bridge elements**
 - Bars
 - Connectors
- Full dentures**
 - Design settings
 - Gingivae
 - Smile Libraries Chart
- Removables**
 - Frames
 - Connectors and Rests
 - Wax profile strips
 - Retentions
 - Stippled waxes
- Digital model**
 - Digital model design
 - Articulator interfaces
- Appliances**
 - Orthodontics Control Panel
 - Positioning Guide

Windows taskbar: 7:02 PM, 1/30/2019, ENG US

Step 4. Click Import Materials.



Dental System Control Panel - Export/Import Materials

Home Back Save Help

Tools > Import/Export

Tasks

- Export materials
- Import materials

Import materials

To import new materials, click on "Import materials" and browse the desired file (extension ".dme"). Then select the materials from the list that appeared on a screen and click "Import". Typically it is recommended to select all. You can also download the materials directly from FTP - go to Site Settings -> Sites page and click "Download Materials".


Import materials

History Of Import Materials

Export materials

Clone system

Step 5. When Windows Explorer is opened, navigate to the file "PREAT_3ShapeLibrary_LabV11.1_ENCRYPTED.dme", select the file and click Open



Step 6. Click Import.



Dental System Control Panel - Export/Import Materials



Tools > Import/Export

Tasks

- Export materials
- Import materials

Import materials

To import new materials, click on "Import materials" and browse the desired file (extension ".dme"). Then select the materials from the list that appeared on a screen and click "Import". Typically it is recommended to select all. You can also download the materials directly from FTP - go to Site Settings -> Sites page and click "Download Materials".

Import materials

History Of Import Materials

Export materials

Clone system

Import and update materials

Please select materials to import

☒ Available materials

☒ Implant systems

☒ PREAT 3i Certain Library (.LAB)

☒ PREAT Ankylos Library (.LAB)

☒ PREAT Astra Library (.LAB)

☒ PREAT Branemark Library (.LAB)

☒ PREAT Camlog Screw-Line Library (.LAB)

☒ PREAT Hiossen Library (.LAB)

☒ PREAT MegaGen AnyRidge Library (.LAB)

☒ PREAT Multi-Unit Library (.LAB)

☒ PREAT Nobel Active/Conical Library (.LAB)

☒ PREAT Nobel Tri-Lobe Library (.LAB)

☒ PREAT Straumann Bone Level Library (.LAB)

☒ PREAT Straumann Tissue Level Library (.LAB)

☒ PREAT Zimmer Screw-Vent Library (.LAB)

☒ Manufacturing processes

☒ Materials

Library name

Provider

Library type

Version

Library creation date


Description

Import


Cancel



Step 7. Click Save. To verify a successful import, click on Implant Systems in the Abutments menu.



Step 8. Verify Preat Implant systems are available in the Implant System window.



Thank You



If you still require technical assistance, please
contact PREAT at:

1-800-232-7732 or technical@preat.com